



ClaraLuna

reference guide



Technical contact :

Manuel MESTRE - Consultant, PhD.
65 avenue de Gournay, 94800 Villejuif, France
+33147260918 / +33608533675
manuel.mestre@free.fr



Commercial contacts :

Didier PASQUELIN - Manager
12 avenue des Prés, BL 128, Montigny le Bretonneux
78059 Saint Quentin en Yvelines, France
tel: +33 1 30 64 05 62
fax: +33 1 30 64 06 38
email: claraluna@optics-concept.com
web: www.optics-concept.com



Peter FLUG - Manager
Fineoptix GmbH
Schmalheck 18, 35625 Hüttenberg, Germany
tel: +49-6441-210 4910
fax: +49-6441-210 4091
email: info@fineoptix.com
web: www.fineoptix.com

Table

First section of this guide : Getting started in one hour

Second section : Advanced documentation

Some features are available only in the "Enhanced version" of ClaraLuna and denoted by a green-circled index: 1

1	GETTING STARTED WITH AN EXISTING FRINGE IMAGE	8
1.1	Getting started – Overview : ClaraLuna inputs and outputs	9
1.2	Getting started – Step 1 : open an image from the disk	10
1.3	Getting started – Step 2 : Edit the Mask	11
1.4	Getting started – Step 2 : Edit the Mask	12
1.5	Getting started – Step 3 : Edit the Options	13
1.6	Getting started – Step 4 : Compute - Analyse results	14
1.7	Getting started – Step 5 : Consider Results files	15
2	GETTING STARTED WITH VIDEO IN EASY-TO- USE MODE	16
2.1	Getting started with Video in Easy Mode – Step 1 : Switch to Easy Mode	17
2.2	Getting started with Video in Easy Mode – Step 2 : Open Video window	18
2.3	Getting started with Video in Easy Mode – Step 3 : Adjust video settings - Grab image	19
2.4	Getting started with Video in Easy Mode – Step 4 : Rename image and Launch	20
3	USING THE ADMINISTRATOR SETUP	21
3.1	Administrator setup : Workspace folder (creating - backing up)	22
3.2	Administrator setup : ISO/DIN options	23
3.3	Administrator setup : Customer data on control reports	24
3.4	Administrator setup : Color scale - Demonstrator setup	25


4	CLARALUNA MAIN WINDOW - MAIN MENU	26
4.1	Main window – Menu : File	27
4.2	Main window – Menu : Setup	28
4.3	Main window – Menu : Analysis	29
4.4	Main window – Menu : Display	30
4.5	Main window – Menu : Help & Tutorial	31
4.6	Main window – Menu : About ClaraLuna	32
	Main window – Status panel & Mode toolbar	33
5	USING THE FILE EXPLORER	34
5.1	ClaraLuna File Explorer: Overview	35
5.2	File Explorer : Shell tree	36
5.3	File Explorer : File list	37
5.4	File Explorer window : Saving a file	38
5.5	File Explorer window : Other controls	39
5.6	File Explorer window : Other controls	40
5.7	File Explorer window : Other controls	41
5.8	ClaraLuna File Explorer vs Windows Explorer : keyboard and mouse commands	42
5.9	ClaraLuna File Explorer vs Windows Explorer : popup menus	43
6	USING THE GENERAL SETUP	44
6.1	General Setup : overview	45
6.2	General Setup - View options	46
6.3	General Setup - Automatic choices (part 1)	46
6.4	General Setup – Automatic choices (part 2)	47
6.5	General Setup – Models of Options and Mask	48
6.6	General Setup - On opening the File Explorer / On creating a new project	49
6.7	General Setup - On opening a new project	50
6.8	General Setup - Buttons	50
6.9	General Setup - Users	51

7	USERS	52
7.1	Login and password access - Users maintenance	53
7.2	User's personal folders and Setup	54
7.3	Access users maintenance from the General Setup	55
8	USING THE MASK EDITOR	56
8.1	Mask Editor - Overview	57
8.2	Mask Editor - Model of Mask	58
8.3	Mask Editor - Building a Mask - Contours	59
8.4	Mask Editor - Drawing contours (ellipse, rectangle)	60
8.5	Mask Editor - Drawing contours (polygons)	61
8.6	Mask Editor - Selecting contours for editing	62
8.7	Mask Editor - Editing contours	63
8.8	Mask Editor - (X,Y) scale	64
8.9	Mask Editor - Master and Disciple Contours ("fiducials") - Enhanced version only	65
8.10	Creating QED 4-Spots Mask in the Mask Editor	66
8.11	Adjusting the QED Mask	67
8.12	Launching and monitoring computation	68
9	CONFIGURING THE PROJECT OPTIONS	69
9.1	Project Options - Overview	70
9.2	Project Options - Model of Options	71
9.3	Project Options - Memos - Reconstruction	72
9.4	Project Options - Optical parameters	73
9.5	Project Options - ISO-DIN specifications	74
9.6	Project Options - Saving Results - Saving / Loading Options	75
9.7	Project Options - Viewing Results	76

10	THE ISO-DIN RESULTS WINDOW	77
10.1	The ISO-DIN results window - Overview	78
10.2	The ISO-DIN results window - 2D - 3D graphs	79
10.3	The ISO-DIN results window - Cross sections	80
10.4	The ISO-DIN results window - Simulated (synthetic) fringes	81
10.5	The ISO-DIN results window - ISO specifications and results	82
10.6	The ISO-DIN results window - Advanced ISO results with confidence intervals	83
10.7	The ISO-DIN results window - DIN specifications and results	84
10.8	The ISO-DIN results window - Zooming	85
10.9	The ISO-DIN results window - Informations and Control buttons	86
11	VIDEO FRAME GRABBING	87
11.1	Video frame grabbing - Entering the Video window	88
11.2	Video frame grabbing - Control buttons	89
11.3	Video frame grabbing - Freeze / Animate video	89
11.4	Video frame grabbing - Computing the displayed image	90
11.5	Video frame grabbing - Saving and processing the displayed image	91
11.6	Video frame grabbing - Saving the displayed image with automatic numbering	92
11.7	Video frame grabbing - Video setup	93
12	EASY PRODUCTION MODE	94
12.1	Easy production mode : Open Mode and create Host	95
12.2	Easy production mode : Open Mode and create Host	96
12.3	Easy production mode : Enter specifications and grab Model image	97
12.4	Easy production mode : Define Series Mask	98
12.5	Easy production mode : Computing - Results window	99
12.6	Easy production mode : See Series results	100
12.7	Easy production mode : Recall previous series - Delete Series	101

13	ADVANCED PRODUCTION MODE	102
13.1	Advanced Production Mode - Purpose	103
13.2	Advanced Production Mode - Entering the Series window	104
13.3	Series of components - Creating a new Series - Folders	105
13.4	Series of components - Creating a new Series - Model	106
13.5	Series of components - Using a Series in Production Mode (1)	107
13.6	Series of components - Using a Series in Production Mode (2)	108
13.7	Series of components - Series Results : seeing statistics, printing, editing.	109
13.8	Series of components - Editing results data	110
13.9	Series of components - Tools	111
13.10	Series of components - Sharing Series between Production Modes	112
14	RADIUS OF CURVATURE	113
14.1	Radius of Curvature for a single part : Enter bench settings	114
14.2	Radius of Curvature for a single part : Enter part settings	115
14.3	Radius of Curvature for a single part : Computing	116
14.4	Radius of Curvature for a single part : Results in ISO/DIN window	117
14.5	Radius of Curvature in Production Mode : Loading series	117
14.6	Radius of Curvature in Production Mode : Computing - Series control report	118


1 GETTING STARTED WITH AN EXISTING FRINGE IMAGE



This chapter gets you started in 15 minutes :

- Open an interferogram from a disk
- Define a Mask and Options
- Compute and read results

1.1 Getting started – Overview : ClaraLuna inputs and outputs



ClaraLuna treats static fringes.

It is a fair alternative to phase shifting devices, when this technique is not relevant for technical or economical reasons.

Advanced Image and Signal Processing provide exceptional performances :


- If the measurement conditions are difficult, ClaraLuna is highly robust to noise.
- If the conditions are good, ClaraLuna can provide Z resolution and precision better than $\lambda/100$, due to linear prediction algorithms : the wavefront reconstruction yields the mathematical expectation, i.e. the most probable shape of the wavefront with respect to the fringe data.

1.2 Getting started – Step 1 : open an image from the disk


Action

Result


1 In Main Menu Setup, click Demo Mode




ClaraLuna File Explorer shows up







2-3 Click Demo01_Defoc.jpg (2) then "Launch" button (3)







Demo01_Defoc is loaded on launching pad







1.3 Getting started – Step 2 : Edit the Mask

Action	Result
<p data-bbox="159 268 739 304">4 To edit the Mask, click button (4)</p> 	<p data-bbox="1169 268 1512 304">The Mask editor shows up.</p> 
<p data-bbox="159 815 593 852">5 Click "Select" tool (5)</p> 	<p data-bbox="1169 815 1579 852">The cursor turns to a square</p> 







1.4 Getting started – Step 2 : Edit the Mask

Action	Result
<p>6 Select contour by clicking it (6)</p> 	<p>The selected contour is shown in finer line with handles</p> 
<p>7 Click and move handles to translate (7), resize (8), rotate (9)</p> <p>...</p> <p>10</p> 	<p>Accept modifications by button (10) and return to Launching pad</p> 

1.5 Getting started – Step 3 : Edit the Options

Action	Result
<p data-bbox="156 279 716 311">11 Now enter the Options window</p> 	<p data-bbox="1164 279 1937 311">The Options window shows up with default parameters</p> 
<p data-bbox="156 790 996 853">12 Adjust a parameter, for instance the measurement wavelength (12) then validate (13)</p> <p data-bbox="156 853 235 885">13</p> 	<p data-bbox="1164 790 1691 821">Now you're back, ready for launching</p> 

1.6 Getting started – Step 4 : Compute - Analyse results

Action	Result																																													
<p>14 Click button "Launch"</p> 	<p>Computation begins : after a few seconds, you get the ISO/DIN Control Report</p>  <table border="1" data-bbox="1534 351 2060 742"> <thead> <tr> <th colspan="2">SPECIFICATIONS ISO 10110 :</th> <th colspan="3">ISO DATA</th> </tr> <tr> <th colspan="5">3/A/[B/C] RMS<D</th> </tr> </thead> <tbody> <tr> <td>PV</td> <td>1.000</td> <td>[0.926</td> <td>0.931</td> <td>0.934]</td> </tr> <tr> <td>A (sagitta error)</td> <td></td> <td></td> <td>0.73</td> <td></td> </tr> <tr> <td>B (irregularity)</td> <td>0.400</td> <td>[0.436</td> <td>0.441</td> <td>0.444]</td> </tr> <tr> <td>C (rotat. irreg.)</td> <td></td> <td></td> <td>0.021</td> <td></td> </tr> <tr> <td>D : RMS t (total)</td> <td></td> <td>[0.220</td> <td>0.220</td> <td>0.220]</td> </tr> <tr> <td>RMS i (irreg.)</td> <td></td> <td></td> <td>0.09</td> <td></td> </tr> <tr> <td>RMS a (asym.)</td> <td></td> <td></td> <td>0.09</td> <td></td> </tr> </tbody> </table>	SPECIFICATIONS ISO 10110 :		ISO DATA			3/A/[B/C] RMS<D					PV	1.000	[0.926	0.931	0.934]	A (sagitta error)			0.73		B (irregularity)	0.400	[0.436	0.441	0.444]	C (rotat. irreg.)			0.021		D : RMS t (total)		[0.220	0.220	0.220]	RMS i (irreg.)			0.09		RMS a (asym.)			0.09	
SPECIFICATIONS ISO 10110 :		ISO DATA																																												
3/A/[B/C] RMS<D																																														
PV	1.000	[0.926	0.931	0.934]																																										
A (sagitta error)			0.73																																											
B (irregularity)	0.400	[0.436	0.441	0.444]																																										
C (rotat. irreg.)			0.021																																											
D : RMS t (total)		[0.220	0.220	0.220]																																										
RMS i (irreg.)			0.09																																											
RMS a (asym.)			0.09																																											
<p>Click "Tilt subtracted" or "Whole with tilt" or "Irregularity" to view different ISO components of the wavefront (15).</p> <p>15 </p> <p>16 </p> <p>Click "3D" or "Simulated fringes" for different graphs (16).</p>	<p>Finally, close window (17).</p>  																																													


1.7 Getting started – Step 5 : Consider Results files

Action


Result

14 Reopen the File Explorer

- either by menu File/Open image from Disk (14.1)
- or by Tool bar icon "Open folder" (14.2)
- or by clicking anywhere on Main window (14.3)



16



Do a mouse right-click on file name for a full screen preview, printing to paper or pdf, sending by e-mail etc...

15


Notice that the Results subfolder has been created in the same Host folder than the interferogram Demo01_Defoc.jpg.

This subfolder has the same name, plus "_Result" :


In the File Explorer, open subfolder Demo01-Defoc_Result (15.1)

See the Result files generated by the previous computation (15.2)

Among others, the Control Report is ready for reading and printing (15.3).





2 GETTING STARTED WITH VIDEO IN EASY-TO- USE MODE






This chapter shows how to :

- **Acquire an image through Video frame grabber**
- **Save the image**
- **Get ready for launching as in previous chapter**

2.1 Getting started with Video in Easy Mode – Step 1 : Switch to Easy Mode

Action	Result
<p>1 In main menu Setup, enter (or stay in) Demo Mode This mode provides a simulation of video framegrabbing even if no video card is installed.</p> <p>2 Then, in main menu Setup, enter Easy-to-use Mode</p>  <p>The screenshot shows a menu with the following items and keyboard shortcuts:</p> <ul style="list-style-type: none"> Administrator Setup (Maj+Ctrl+A) General Setup (Maj+Ctrl+G) Options of Current Project (Maj+Ctrl+O) Printer Setup (Maj+Ctrl+I) Regular Mode (F2) Easy-to-use Mode (F3) - circled in blue One Shot Mode (F4) Easy Production Mode (F5) Advanced Production Mode (F6) QED Production Mode (F7) QED Spots Mode (F8) Video Sequence Mode (F9) Gratings Mode (F10) Demo Mode (F11) - circled in blue 	<p>Now you are in Demo+Easy Mode</p>  <p>The Easy Mode provides the following behaviour :</p> <p>OPEN FILE EXPLORER: <input checked="" type="radio"/> in last opened folder</p> <p>When the File Explorer opens, the folder is the last used.</p> <p>MASK <input checked="" type="radio"/> Last used Mask</p> <p>When an interferogram is launched, the initially prompted Mask is the last used</p> <p>OPTIONS <input checked="" type="radio"/> Last used Options</p> <p>and the initially prompted Options are the last used.</p>

2.2 Getting started with Video in Easy Mode – Step 2 : Open Video window

Action	Result
<p>3 In main menu File, enter Acquire Video Image</p> 	<p>Now you are in the Video window</p> <p>Because of the Demo Mode you selected, the image does not actually come from a video input. However, the behaviours of most controls are identical to those of real video, without demo.</p>
<p>Alternatively you can call Video from the File Explorer</p> <p>First open File Explorer by menu (3.1) or toolbar button, or click on main window. When in File Explorer, click button (3.2)</p> 	



2.3 Getting started with Video in Easy Mode – Step 3 : Adjust video settings - Grab image

Action

Result

4

In the Video window, click "Tools" button (4)



In real video, you can adjust the video standard of your camera, and image Brightness/Contrast.


"Saturation on" helps you detect whether some areas on the image are too dark or too bright.

(Because of the Demo Mode, these controls are deactivated in this tutorial)

Click "Ok" button to close video tools.


5

Then click "Save" button (5)







The File Explorer opens in the last opened folder (6)


The grabbed image is visible (7). It is called TempVideo.jpg (8)



2.4 Getting started with Video in Easy Mode – Step 4 : Rename image and Launch

Action	Result
<p>9 ... 12</p> <p>Click "Save as" edit box and enter a new file name, for instance "Test1", then press "Enter" key or click "Save" button (10).</p> 	<p>Now the grabbed video image is called "Test1.jpg" (11) ready for launching (12)</p> 
<p>12</p> 	<p>MASK</p> <p><input checked="" type="radio"/> Last used Mask</p> <p>OPTIONS</p> <p><input checked="" type="radio"/> Last used Options</p>  <p>From now on, proceed as in previous chapter "Getting started with an existing fringe image".</p> <p>Recall that, due to the Easy Mode, the initially prompted Mask and Options are the last used ones.</p> <p>If necessary, modify them (13,14), then launch computation (15).</p>

3 USING THE ADMINISTRATOR SETUP



In this chapter, set the general behaviour of ClaraLuna common to all users:

- Define folder ClaraLuna Workspace
- Set important ISO/DIN output options
- For showing on control reports, enter color scale, customer logo and info

3.1 Administrator setup : Workspace folder (creating - backing up)

1 Enter Administrator setup
by Main menu Setup> Administrator setup
or shortcut Shift+Ctrl+A

2 Select Customer Setup
in the outer tabsheet




3 Select “ Workspace folder ” in the inner tabsheet

- 4 • The default location for “ClaraLuna Workspace” folder is on the usual data partition D:\, if any, and on C:\ if not.
- 5 • If the current “Workspace folder” is already defined in ClaraLuna Setup file, and/or it can be found at the default location on the disk, you will prompted to confirm it (5).
- If “ClaraLuna Workspace” folder already exists, it is recommended not to create a new one, in order to avoid confusion or loss of data and setup information.
- 6 • If the “Workspace folder” is not defined in ClaraLuna Setup file, and cannot be found at the default location on the disk, you will be asked to create it at the prompted default location, or at a custom location - to be selected on the File Explorer (6) in the same window -

The best Backup policy is to use the Workspace folder for all your projects (saving and computing fringe images).

7 For backing up all you data (i.e. all configuration files, computed projects, masks...) just zip up the whole “ClaraLuna Workspace” folder : use any explorer - Windows’ or ClaraLuna’s - right-click on the folder icon and call any file archiver such as Winzip, or the freeware IZarc distributed with ClaraLuna, and available for installation in folder C:\Program Files \ClaraLuna\Utilities

Save this archive. To restore, just unzip the archive in the drive where the Workspace is, or should be. If you choose unzip option “overwrite existing files”, you will get back all your previous settings and data, without deleting newer files.



3.2 Administrator setup : ISO/DIN options

1 Select sheet (1)

2 Making results more specification-sensitive

As described in Section 10.6, the tolerance (if any) of a given ISO/DIN parameter is shown on the ISO/DIN control report in green or red according to the result. Orange is for parameters with a confidence interval (3).

The setting (2) shifts by s% (a security coefficient) the values that control the specification colors :

- the measured values and the tolerances are unchanged
- the thresholds for color orange or red increase
- as a result, in case no critical parts and tolerances, the user is warned if the risk of being out of tolerance increases


4 In contrast, the setting (4) “Making ISO/DIN results smaller” allows to cut a given percentage (e%) of the surface extreme Z data, for computing ISO/DIN parameters.

This makes Peak to Valley parameters (PV, B and C) much more robust with respect to noise and random fluctuations. Such an operation makes sense from an optical point of view, since the optical properties of components usually rely on average values (RMS, Strehl ratio...), and seldom on extreme values (such as maximum local slope for power lasers).

This setting truncates the histogram of Z values by e/2 % at both ends, then compute Peak to Valley and RMS on the remaining data.

- ### 5 CAUTION : This setting should be done only by qualified users, aware both of :
- the technical-commercial meaning of the part drawing,
 - and how these choices may impact the reception control by the final customer.




USING THIS SETTING SHOULD BE DONE ONLY AFTER VALIDATION BY THE FINAL CUSTOMER, and the percentage allowed for cutting off should be written on the part drawing.



SPECIFICATIONS ISO 10110 :		ISO DATA		
3/A/(B/C) RMS _{x<D}				
PV	60	[41	43	53]
A (sagitta error)		1.4		
B (irregularity)	50	[40	43	53]
C (rotat. irreg.)		3.0		
D : RMS t (total)	6.0	[6.0	6.1	7.6]
RMS i (irreg.)		5.9		
RMS a (asym.)		6.0		


3.3 Administrator setup : Customer data on control reports

- 1 **Select sheet (1)**
- 2 **Open bitmap (.bmp) file of firm logo.**
For better graphic results, the image size should be : width \leq 180 and height \leq 80, and preferably as large as possible. If the original image is larger, it will be reduced to fit. As a result the fine image details, especially characters, will be blurred.
- 3 The opened image will show in (3)
- 4 **Write the firm info in the edit box (4)**
Name, address, phone, web...
- 5 **Once the firm logo and info are set** you may want to show them or not on the control reports, by selecting in (5) If “Show “ is selected, logo and info will be shown on the control reports.
- 6-8 **The control reports are created when:**
 - Printing (6) the ISO/DIN form, in the present stage, to paper or pdf.
 - Saving (7) the ISO/DIN form, in the present stage, to jpeg.
 - Closing (8) the ISO/DIN form : it will be saved, in its last stage, to jpeg.
- 9 **The control reports will show the logo and info in (9)**




3.4 Administrator setup : Color scale - Demonstrator setup

- 1 Select sheet (1)
- 2-3 Select color scale
- 4 Options for demonstrators can be set in tabsheet Demonstrator Setup
- 5 • Show/hide Demo icon on Modes toolbar at the right hand side of main window
- 6 • On calling the Demo Mode, call/do not call File Explorer
- 7 • As default setting, all changes in General Setup done while Demo Mode is on, are discarded when exiting Demo Mode. For demonstrators, it may however be useful to keep settings, such as Series.
- 8 • Diffraction gratings is a special application, whose demo is time consuming. Keep unchecked in most cases, unless you really need to demonstrate Gratings Mode.
- 9 Developer Setup is only for maintenance, and is password protected.



4 CLARALUNA MAIN WINDOW - MAIN MENU




This chapter tackles :

- **Main window**
- **Overview of menus and toolbar**

4.1 Main window – Menu : File

- 1 Exit**
Closes ClaraLuna.
(In latin, "exit" = "he goes out")

- 2 Acquire video image**
Opens Matrox frame grabbing window
- 3 Acquire Twain image**
Opens frame grabbing window for all Twain compliant devices : camera, scanner...
- 4 Open ClaraLuna File Explorer**
This is the main tool for opening image files from a drive (hard disk, CD-DVD, USB device, local network, or the Web).
Save, copy, rename, change file format.
Define file automatic numbering or dating, with custom names.
For comfortable navigation, all operations are memorised : visited folders, opened images, computed projects...
- 5 Shortcuts to 10 previous projects**
- 6 Keyboard shortcuts for commands**
The usual Windows shortcuts are used (Alt F4 to exit, Ctrl F4 to close active window), but most ClaraLuna's shortcuts are Ctrl+Alt+ some letter).
- 7 Toolbar**
Click images for the most usual commands.



4.2 Main window – Menu : Setup

1

General setup

Here you will configure the general behaviour of the software, effective on all projects.

2

Options of current project

Each interferogram, when computed, becomes a Project with a specific configuration, optical parameters, and specifications. The Project Options window is the project's headquarters.

3

Administrator options

Used only for special applications and maintenance.

4

Printer Setup

For sending the Results windows (ISO, Zernike, slopes...) to a printer for paper copy, or PDF file.

5

Regular Mode

Normal way of working with ClaraLuna. Equivalent to driving a car with a manual gearbox.

6

Easy-to-Use Mode

Could be of help for a start : No configuration needed, everything (mask, File Explorer) stays exactly in the state it had been left.

7

One Shot Mode

Calls video and computes without saving, for instant part testing.

8

Easy Production Mode

Minimum configuration needed. Straightforward and fast for large industrial series and low-skilled users

9

Advanced Production Mode

Assistant for advanced configuration when using ClaraLuna for a series of parts. Just like driving a car with an automatic gearbox.

10

QED Production Mode and Spots Mode

For special applications with QED -MRF[®] machines

11


Video Sequence Mode

Programmable time series of video acquisition : N batches of n images are grabbed, with custom delays between batches and between images, followed by computing.

12

Gratings Mode

For special applications to interferometric measurement of diffraction gratings.



13

Demo Mode

Emulates a video fringe frame grabbing. Everything else is actual computation. A panel of preset demos is prompted, with typical interferograms ranging from exceptionally high quality and precision, to rather bad conditions (noise, low contrast...)

4.3 Main window – Menu : Analysis

- 1 Zernike Polynomials**

Access the Zernike window after computing an interferogram : interactive analysis, removal of any set in the 36 first Z.
- 2 Legendre Polynomials**


Advanced feature, suitable for rectangular apertures ; fitting only enhanced versions of ClaraLuna.
- 3 Wavefront slopes analysis**

Specifications compared to ptv and rms of slope maps (enhanced versions).
- 4 PSF - MTF – Strehl ratio**

In enhanced versions, on request.
- 5 (Re)compute current project**


Launches computation for the project whose window is focused. Recomputes with new options if the previous result is not satisfactory.
- 6 (Re)compute all open projects**

Same effect on each project window, already computed or not.





4.4 Main window – Menu : Display

- 1 Cascade**
Displays the open project windows
- 2 Tile**
Opens the open project windows
- 3 Close active window**
Closes the active project window
- 4 Close All**
Closes all open projects windows.
- 5-6 Previous / Next window**
Selects and brings to front the next open projects window
- 7 Minimise all windows**
Reduces the open projects windows
- 8 Expand all windows**
Restores the size of the open projects windows
- 9 Show tool bar**
Shows/hides toolbar at the bottom of the window
- 10 Language**
English - French available.
Chinese is on the way.
Others on request



4.5 Main window – Menu : Help & Tutorial

- 1 **This Reference guide**
“Getting started” is the first 2 chapters
- 2 **Help on / help off**
Integrated help (in English or French only) :
interactive popup windows provide detailed
explanations about what's pointed by the mouse.
- 3 **Optimising the interferogram**
What fringes are from a geometrical point of view.
How ClaraLuna reconstructs the part surface from
the fringe data. How to get poor data and results,
and how to get better ones.
- 4 **More information on special functions**
A Shockwave presentation to ClaraLuna's main
features
- 5 **Calling Demo mode**
Same Mode than through menu Setup>Demo
Mode, with an additional “how to use” message.
- 6 **Flash movie tutorials**
- 7 **Scientific background** for ClaraLuna's
wavefront reconstruction algorithms
- 8 **Testbed** for calibrating ClaraLuna with Zygo
reference Phase-shifting device and software,
giving evidence of $\lambda/100$ precision and
resolution.
- 9 **Paper** by J. Briers in Journal of Optics, reporting a
inter-laboratory comparison : the same sample of
optical components have been sent to various
labs. The range of the results is astonishing. This
helps putting in perspective the precision in
optical metrology...



4.6 Main window – Menu : About ClaraLuna

1 Informations on this version of ClaraLuna

List of all the functions available in the full version of ClaraLuna, and which ones are licensed to you.

The license is included in the Sentinel Key provided with the software. All executables of a given version (4.38 etc) are identical, regardless of the license.

For upgrading the license and purchase additional functions, click :

2 Request for upgrading licensed features

You will need to send a “.req” file, receive a “.upw” file, and play it through the Sentinel utility provided.

3 Contact


Commercial and technical contacts : asking for help or tutorial, reporting (dis)satisfaction, upgrading towards enhanced functions...

4 License agreement

Legal agreement between the author and the user. Though it is not fun, please read it ! The main idea is that you cannot charge the author or his agents for any trouble or consequence due to using this software ...

5 Credits

Who was involved in the creation of this software, as well as used libraries and utilities.



4.7 Main window – Status panel & Mode toolbar

1 Status panel

At the right hand side of the main window, the status panel shows the current User (2) (see Section 6.9 and Section 7) and the current Mode (3).

2

3

4 Mode toolbar

The icons (4) are shortcuts to the Modes listed in the Setup menu (5), described in section 4.2.


5

6 If the selected Mode is one of the three Production Modes (Easy, Advanced or QED), the current Series of parts is shown on the status panel (6). The keyboard shortcuts for the Modes are the F keys (5): F2...F9 (see Setup menu).

The diagram illustrates the main window's status panel and mode toolbar. It shows a screenshot of the software interface with numbered callouts (1-6) pointing to specific elements. The status panel (1) displays the current User (2) as 'Paul' and the current Mode (3) as 'Regular'. The mode toolbar (4) contains icons for various modes, including 'Regular Mode', 'Easy-to-use Mode', 'One Shot Mode', 'Easy Production Mode', 'Advanced Production Mode', 'QED Production Mode', 'QED Spots Mode', 'Video Sequence Mode', 'Gratings Mode', and 'Demo Mode'. A setup menu (5) is shown below the toolbar, listing keyboard shortcuts for these modes. The status panel also displays the current Series (6) as 'DemoSeries2' when in a production mode.

General Setup	Maj+Ctrl+G
Options of Current Project	Maj+Ctrl+O
Administrator Setup	Maj+Ctrl+A
Printer Setup	Maj+Ctrl+I
Regular Mode	F2
Easy-to-use Mode	F3
One Shot Mode	F4
Easy Production Mode	F5
Advanced Production Mode	F6
QED Production Mode	F7
QED Spots Mode	F8
Video Sequence Mode	F9
Gratings Mode	F10
Demo Mode	F11

5 USING THE FILE EXPLORER



In this chapter :

- Open image files from a drive (hard disk, CD-DVD, USB device, local network, or the Web).
- Save, copy, rename, change file format.
- Define file automatic numbering or dating, with custom names.

5.1 ClaraLuna File Explorer: Overview

1

Shell tree : Drives and folders

Displays a view of the whole PC working station

2

Shell list : Files

Displays the content of the open folder: subfolders and files.

The right panel shows the selected file if it is a jpg, bmp, int or dat image.

3

"Save as" :


Copies the visible image under a new name.

4

Buttons

ClaraLuna File Explorer vs Windows Explorer :

The ClaraLuna File Explorer is Windows Explorer compliant. The behaviour and keyboard commands are identical. You can use drag and drop within the File Explorer, and also between the File Explorer and Windows Explorer. Right clicking in the File Explorer makes it possible to call the Windows Explorer function "search", or open a file with the associated software....



5.2 File Explorer : Shell tree

1-2

Select folder

Click a folder to select and open it. You can as well access a folder by typing or pasting its path in the combo (3) and pressing "Enter" key. This function is useful in particular when opening a file with a very long and known path. This happens after a Search action in the Windows Explorer, by copying the found file complete path from the Windows Explorer (2) to the Clipboard.

3

History of open folders

Click list and select a previously visited folder.

Selecting a folder in the list opens it in the shell tree (1) and puts its name on top of list (3). Selecting a new folder in the shell tree (1) also puts its name on top of list (3).

4

Browse folders history

When selecting a new folder in shell tree (1), the list of the previously visited folders is fixed. In the example, you visited the subfolder \1, then \2, then \3. Clicking the left arrow (3) takes you back in the history list : \2, \1, then \3 again. The right arrow will visit these folders the other way round.

As long as you do not add a new folder, even if you modify the order of the list (3) by selecting an item, the arrows (4) will follow the effective time order of the folder opening.

When closing the window, the history list is saved.

5

Headlines of PC working station

Select one of the Desktop folders or one of the drives (Hard disks, CD-DVD, removable storage media)

6

Create subfolder

Click to create new subfolder in the selected folder, then rename it.

7



Parent folder

Click to select the "parent" folder of the current open folder (on the image, \ClaraLuna Demo\ is the parent folder of \1). Keyboard command: press key "delete left of cursor" (left arrow on top right of keyboard).

8

Delete selected folder

... and all its content. Keyboard command: press key "delete" (Del).



5.3 File Explorer : File list

1

History of displayed images


All images shown in the viewer have their path memorised in this list. Click (1) and select an item to open its folder and view the image.

2

History of projects

When launching a fringe image, a Project is created with the image name: see chapter "Computing an interferogram".

Button (2) opens the history of previously created Projects.



3-7


Open a file by entering its name

When clicking combo (3), combo (4) and buttons (5) (6) appear. Enter a file name ("Demo 15") or a full path ("C:\ClaraLunaData\Demo15.jpg").

If you click a file in the list (7), its name is copied in combo (3) for editing.

If you enter just a file name ("Demo 15"), the extension (".jpg") will be added according to selection in combo (4). Otherwise the extension you write ("Demo 15.bmp") will be used.

Finally, click (5) to open file, or (6) to cancel .



5.4 File Explorer window : Saving a file

1-3

Saving the displayed image as ...

The image shown in the viewer can be saved under a new name and/or a new format. Click combo (1) and enter new name. Check or modify format (2) Click button (3).

Automatic file naming and numbering

4

Select in the list a format for file names.

5-6

Choosing "New..." or "...Comment..." will let you enter some text, for instance "Mirror".

7

Validate on green tick.

8

Your custom file name is included in the list. It is saved and can be reused later.

9

Click button to save the displayed image.


10

The first image is saved under "Mirror000.jpg".
The second under "Mirror001.jpg".

11

Click red cross to delete custom file name.

This function is very comfortable when acquiring series of data from a video frame grabber.



5.5 File Explorer window : Other controls

1

Showing legends in the File Explorer

Check box to show the legends.

2

Even if the legends are hidden, most controls are described by a "hint", ie a short comment popping up when then mouse cursor stops over the control during one second.

Central buttons

3

Launches computation : the displayed image will now be turned to a Project, with a Results folder created in the image folder, containing all the related files. But what you will see in the first place is the File Explorer being closed, then the interferogram appearing in the Project window.

You get the same effect by clicking the fringes image in the viewer, or double-clicking the image name in the file list.

4

Opens Video window (same effect than opening it from the Main window).

5





Closes File Explorer.

6

Opens Help.

Save settings

When clicking one of the exit buttons - launch (3), video (4) or cancel (5) - all the File Explorer settings will be saved before closing the File Explorer, and stored into the ClaraConfig.ini file.






5.6 File Explorer window : Other controls

1 Show thumbnails

2 Show folders

3 Hide folders



5.7 File Explorer window : Other controls

1 Current user's personal folder

2 Click (1) for reaching the current user's personal folder (2) located in ClaraLuna Workspace\Users


3 This home folder contains subfolders :

- Config (holding the user's setup ".ini" files)
- Masks (holding a copy of all the Masks created and/or used by the current user)
- Options (holding a copy of all the Project Options files created and/or used by the current user)
- Projects (if the user wants to use this folder instead of the common projects folder ClaraLuna Workspace\Users)

Edit Export folder for .dat files

This folder, typically located on the local network, contains the files with format Zygo .dat sent to the QED machine - or for other export use.

- 4
- 5
- 6
- Open dialog
 - Choose action
 - Accept and close dialog



5.8 ClaraLuna File Explorer vs Windows Explorer : keyboard and mouse commands

1

ClaraLuna File Explorer is Windows Explorer compliant

... as mentioned before. The behaviour and keyboard commands are identical.

2

Keyboard commands

based on the standard windows keyboard shortcut

3

Mouse commands

Use Ctrl + drag and drop within the File Explorer, and also between the File Explorer and Windows Explorer


Action

Delete selection
Copy selection
Cut
Paste
Select all
Go to parent folder

2


Keyboard commands

Del
Ctrl + C
Ctrl + X
Ctrl + V
Ctrl + A
Back space



*Move or Copy by Drag and drop
or Ctrl X / Ctrl C / Ctrl V
to/from ClaraLuna File Explorer
from/to Windows Explorer*

2-3



5.9 ClaraLuna File Explorer vs Windows Explorer : popup menus

Mouse commands- Calling Windows Explorer's functions from Claraluna's File Explorer

1

A mouse right click on a folder opens a popup menu for folders.

Selecting "Search" will call the Explorer file search window.

Use whatever programs are installed in the menu : typically an antivirus, or a file compression utility.



2

A mouse right click on a file opens a popup menu for files.

Double-clicking a file name in Claraluna File Explorer launches the computation of an interferogram, provided that the image format is read (jpg,bmp,int,dat) and the image is not a postcard from Paris.

No other action by double-clicking a file is allowed directly from Claraluna File Explorer. However, opening any file associated to a program (say Report.doc) is straightforward : right-click Report.doc, select Open in the popup menu : the file will be opened in MS Word.

Finally, the commands "cut", "copy", "paste" are available as well in these popup menus.



6 USING THE GENERAL SETUP



In this chapter :


- Set the general behaviour of ClaraLuna specific to the current user
- Save, retrieve your personal configuration

6.1 General Setup : overview


Purpose


The General Setup manages ClaraLuna's behaviour, common to all Projects. Please read carefully this section. The Setup window contains :

- 1 **View options**
Some appearance options : toolbar, language...
- 2 **Automatic choices**
On selecting several images
Timeout for launching computation
Single/multiple results folder.
- 3 **On opening the File Explorer**
Choose which folder the File Explorer opens
- On opening a new Project**
- 4 Choose which Mask model is prompted
- 5 Choose which Options model is prompted
- 6 **Options and Mask Models**
Choose the properties Write enabled or read only for the Models
- 7 **Post it**
Write down any note or general message for the user of this configuration.
- 8 **Users**
Define Users, each having personal setup, saved to a personal ini file.
- 9 **Buttons**
for saving the displayed config to a personal file ;
for retrieving your personal setup from that file.



6.2 General Setup - View options

- 1 **Maximize ClaraLuna** : If checked, the software spreads over the whole screen.
- 2 **Show logo** On opening ClaraLuna, the logo is shown then hidden after timeout (2).
- Show toolbar** in main window :
- 3 
- 4 **Show hint**
If checked, most controls will be described by a "hint", ie a short comment popping up when then mouse cursor stops over the control during one second.
- 5 **Tutorial messages**
If checked, some messages pop up : "Caution: if you do this, that will happen..."
Uncheck for regular use, as soon as you get familiar with the software.
- 6 **Customise some commands** : when clicking main form, do nothing / call File Explorer... and similar behaviour options.
- 7 **Language**
French/English presently available. Other languages on request... if possible.




6.3 General Setup - Automatic choices (part 1)

On selecting several images

In the File Explorer, you can select several image files (Ctrl+click, or draw a frame with the mouse over the file names). When launching computation, the files will be computed in a row as :

- Independant projects, or :
- (Enhanced version only) Independant projects followed by an averaging. In this case, before computing, you will be prompted a name for the average. All these Projects results (the image files plus the average) will be stored in separate Results folders created in the same folder as the images, with the related names.

- 7 After computation, checking box(8) will close the Results windows of all the Projects, but not the final average (if any) Result window.



6.4 General Setup – Automatic choices (part 2)


Launch computation automatically

5 When launching a Project computation from the File Explorer, the first step results in showing Project window (5) for validating the Mask and Project Options.

6-8 If box (7) is unchecked, window (5) will be on standby, waiting for you to click any of the buttons.

By checking box (7) and setting a timeout (8) in seconds, the countdown (6) will automatically trigger the computation in the Project window (5).

Note that this automatic option might stay checked, in case of an unexpected interruption during the computation of a series of images (independent or averaged). In this case, just uncheck box (7).




Automatic indexing of Results folders

A fringe image (say "Image.jpg") will be turned into a Project by starting the computation from the File Explorer or the Video window. This means that the folder "Image_Result" will be created beside Image.jpg, containing all the files related to the Project : Mask, Options, ISO results...

You may want to recompute "Image.jpg" with different settings :

9 If option (9) is on, the new Results folder will overwrite the previous one.

10 If option (10) is on, the Results folders will not be overwritten, and will be named with the common name "Image_Result" followed by an incremented number.



6.5 General Setup – Models of Options and Mask

1 Mask and Options Models of a new Project :

Beside the interferogram (i.e. the fringes image itself), the main ingredients of a Project are a Mask and Options. The mask is a graphic object that defines the data/no data areas on the interferogram, as well as the "horizontal" XY scale.

The Options include all the parameters specific to this particular Project (optical data such as wavelenghtes, ISO/DIN specifications, results display and saving).

On opening a new project from the File Explorer or from the Video acquisition window, a Mask and Options will be prompted to the user for validation and/or modification before the actual computation.

These prompted Mask and Options are called "Models", and are similar to document models in word processors.

When launching the computation the validated Mask and Options are integrated in the computation and are saved as files with the name of the Project (= the name of the interferogram) + the extension ".msk" for the Mask and ".opt" for the Options. These files are copied to the Results sub-folder, in the same host folder than the interferogram .

1.1 The Mask and Options Models are retrieved from existing files, that can be chosen by the user : (1.1) for Mask, (1.2) for Options - see details in next Section.


1.2 The Mask and Options Models are independant, and the way of choosing them is independant, too.

If the Mask and/or the Options are modified by the user before computation, the changes can be saved to their respective Model or not, depending on the user's will :

2 Saving to Mask and Options Models :

- 2.1 • Save changes to the Models (set a Write-enabled property)
- 2.2 • Do not modify these models (set a Read-only property).

Note that the latter choice holds for **both** Mask and Options Models.




6.6 General Setup - On opening the File Explorer / On creating a new project

Open File Explorer

When opening it, the File Explorer will display a given folder. Choose :

- 1 • In the last folder opened (you will find it in the same place where you last left it). The displayed image will be the first one in the folder whose name is read by ClaraLuna.
- 2-3 • In a specified folder : browse by button (4) for calling the File Explorer and selecting a folder. This option is the best choice when working with a series of parts.
- 4 • With a previous project : come back to the folder of the last Project that was computed, and display its image.




On opening a new Project : Mask

When launching an interferogram from the File Explorer, a Mask will be created from of Model of Mask, and prompted for validation or modification, before the actual computation. Choose :

- 5 • **No mask** if you want to start from a blank. But you cannot compute without defining at least one aperture.
- 6 • **Automatic selection** : finds an ellipse or a rectangle with axes horizontal and vertical. It works well only for a bright fringe image on a dark background.

All Projects (for instance "Demo1_Defoc.jpg") create a mask file with the same name and extension ".msk", and makes two copies : one in the Project Results folder ("Demo1_Defoc_Result"), another in the common folder "Masks" located in the ClaraLuna program folder. If computing several Projects with the same interferogram, the last mask will overwrite the previous one in the "Masks" folder, but not in the Results folder (provided that you choose the option "Automatic indexing of Results folders") :

- 7 • **NameOfImage.msk read in the "Masks" folder**
- 8 • **NameOfImage.msk in the folder of image**
- 9 • **First file .msk read in folder of image** : drop any msk file beside the image
- 10-11 • **Model of Mask specified by user** : browse and select with button (11). This option is the best choice when working with a series of parts.
- 12 • **A previous mask (by default: the last one)**. Easy to use for a first run on ClaraLuna, but not suitable for working with a series of parts.



6.7 General Setup - On opening a new project

On opening a new Project : Options


The choices for Project Options are very similar to those for Masks.

- 1 • **ClaraLuna default options** calls the built-in factory settings, ie the most common and easy-to-use.
- 2 • **NameOfImage.opt read in the "Options" folder** (see page above, # 7)
- 3 • **NameOfImage.opt in the folder of image** (see page above, # 8).
- 4 • **First file .opt read in folder of image** : drop any msk file beside the image
- 5 • **Model of Options specified by user** browse and select with button (7). This option is the best choice when working with a series of parts.
- 6 • **Options of previous Project**

Options and Mask Models :

- **Save changes to the models:** if checked, the Options Model is write-enabled
- **Do not modify the models : read only**

The same checkbox controls the "read only" behaviour of both the Mask and the Options Models



6.8 General Setup - Buttons

- 1 **Load settings**
Reads a .ini file and loads the settings to the General Setup window.
- 2 **Load default settings**
Loads the built-in factory settings
- 3 **Save settings as**
Creates a .ini file with the settings of the General Setup window.
- 4 **Accept**
Saves the settings of the General Setup window to the default file ClaraConfig.ini, located in the "Options" folder - When starting the program ClaraLuna, this configuration file ClaraConfig.ini will be loaded. Note that that the last config will be loaded for the next session. For changing to another saved config, use button (1).
- 5 **Cancel**
Closes window, discards all modifications done to the settings since opening the General Setup window.
- 6 **Help**



6.9 General Setup - Users

1-2 The current user logs in awhen ClaraLuna opens.

3 The current user's name is shown in edit (1).
Click button (2) to open the Users Maintenance window.

4-5 In the Users maintenance window,

- 3**
 - create new user
 - switch to next user
 - delete user
- 5**
 - change password

Back to general Setup window


- 6**
 - If next user has logged in, the change will be taken into account only if closing the General Setup window by button (4),

6-7 In the Main window status panel the new user will be visible (5)

The current User's name is shown in the status panel, top right of main window

8

9



9 →

